



Yes, you can.™

The Invacare® Top End® Eliminator™ OSR is designed for speed. The super-stiff roll-formed horizontal maintube, combined with a relaxed headtube angle for exceptional stability and aerodynamics and an oversized raked fork for maximum traction, make it the choice of champion needs, like Franz Nielispach for four years straight years in his Eliminator OSR, Nielispach has won the most prestigious wheelchair race of all, the Boston Marathon. So if you have your sights set on winning, the Eliminator OSR will help deliver your best, just like it does for Franz.



The Invacare® Top End® Eliminator™ OSR featuring a relaxed head tube angle for a exceptional stability and front wheel cornering.



Oversized rake for maximum speed and aerodynamic steering for optimum control.



Fixed footrest option is adjustable fore / aft, height and angle.



The precision, track compensator system ensures maximum speed on the road and around each turn of the track.



The Invacare® Top End® Eliminator™ OSR with optional front and rear carbon fiber wheels.

Standard Features

- 6061 T6 oval tube aluminum frame
- Integral wrap-around fenders
- 3 cage designs: V,V-open or I-cage
- Kneeling, amputee or traditional upholstery
- 700C or 26" high-performance wheels
- Aluminum side panels
- Straight axle tube with rear wheel alignment
- Aluminum fork
- Precision track compensator system
- 20" front wheel with tubular tire
- Taped 3/8" diameter aluminum pushrims
- Aero bars
- Powder coat paint

Options Available

- Carbon fiber wheels, front and rear
- Fixed footrest
- Cushion
- Tire coated pushrims
- Cordless speedometer
- Helmet
- 18" front wheel
- Custom paint

Invacare Corporation
www.invacare.com

USA

One Invacare Way
Elyria, Ohio
44035-4190
(800) 333-6900

Invacare Top End Sports and Recreation Products

4501 63rd Circle North
Pinellas Park, Florida
33781-5914
(727) 522-8677
(800) 532-8677
Fax: (727) 522-1007

Canada

5970 Chedworth Way
Mississauga, Ontario
L5R 3T9 Canada
(800) 668-5324

Invacare and Top End are registered trademarks and Eliminator and Yes, you can., are trademarks of Invacare Corporation. Specifications are subject to change without notification.

Top End®
Eliminator™
OSR Racing
Chair

Invacare® Top End® Eliminator™ OSR Racing Chair



Ernst Van Dyk, Team Invacare



Yes, you can.™