



Lift Recliner Specifications

Lift Recliner Specifications



TRADITIONAL SERIES

	Williamsburg	Ambassador	Pioneer	Pioneer
Model Number	PR-747	PR-760	PR-640	PR-641
Size	Medium	Medium	Medium	Medium
Weight Capacity	375 LB	375 LB	375 LB	375 LB
Recline Positions	3	3	3	3
Standard Fabrics	7	7	3	3
Fabric Options	36	36	36	36
Chaise Seat	YES	YES	YES	YES
Back Type	BUTTON	3 PILLOW	SEAM	SEAM
Head Rest Cover	YES	YES	YES	YES
Arm Covers	YES	NO	NO	NO
Battery Back Up	YES	YES	YES	YES
Overall Width	32.5"	31.5"	29.5"	29.5"
Overall Height	44"	42"	41"	41"
Width Between Arms	21"	21"	21.5"	21.5"
Floor to Top of Seat	19.5"	19.5"	19"	19"
Seat Depth	21.5"	21"	21"	21"
Seat to Top of Back	28"	27"	28"	28"
Seat to Top of Arm	6"	7"	8"	8"
Floor to Back Edge of Seat Raised	28"	27"	27.5"	27.5"
Floor to Front Edge of Seat Raised	24"	23"	23"	23"
Seat Angle Raised	15 DEG.	22 DEG.	20 DEG.	20 DEG.
Floor to Top of Foot Rest Reclined	23"	23"	23"	23"
Floor to Top of Back Reclined	36"	34"	34"	34"
Back Angle Closed	100 DEG.	104 DEG.	115 DEG.	115 DEG.
Back Angle Reclined	130 DEG.	140 DEG.	140 DEG.	140 DEG.
Foot Rest Angle Closed	90 DEG.	90 DEG.	90 DEG.	90 DEG.
Foot Rest Angle Opened	180 DEG.	180 DEG.	180 DEG.	180 DEG.
Distance Required from Wall Reclined	18"	19"	18"	18"
Extension reclined from Back of Seat to end of Foot Rest	42"	40.5"	40.5"	40.5"
Life Time Warranty, Chair Frame, Lift Frame, Recline Mechanism	YES	YES	YES	YES
Warranty on Electrical parts	3 YRS.	3 YRS.	3 YRS.	3 YRS.
Pro Rated Years 4 Thru 7	YES	YES	YES	YES
Warranty on Mechanical Labor	3 YRS.	3 YRS.	3 YRS.	3 YRS.
Warranty on Electrical Labor	1 YR.	1 YR.	1 YR.	1 YR.
Removable Back	YES	YES	YES	YES
Customizable *	NO	NO	NO	NO
Left Hand Control Available	YES	YES	YES	YES
Extra Pocket Available	YES	YES	YES	YES
Smart Tek System	YES	YES	YES	YES
Heat Massage Option	YES	YES	YES	YES
Foot Rest Extension Option	YES	YES	YES	YES
MaxiComfort Option	NO	NO	NO	NO



SIGNATURE SERIES

	New Yorker	Regal	Royal	Relaxer	Space Saver
Model Number	PR-759	PR-751TY	PR-752	PR-756MC	PR-906
Size	Medium	Medium	Medium	Medium	Medium
Weight Capacity	375 LB	375 LB	375 LB	375 LB	250 LB
Recline Positions	3	3	3	MaxiComfort	3
Standard Fabrics	8	8	8	8	3
Fabric Options	36	36	36	36	36
Chaise Seat	YES	YES	YES	YES	YES
Back Type	4 PILLOW	3 PILLOW	3 PILLOW	4 PILLOW	SEAM
Head Rest Cover	YES	YES	YES	YES	NO
Arm Covers	NO	NO	YES	NO	NO
Battery Back Up	YES	YES	YES	YES	YES
Overall Width	37"	36"	32.5"	36.5"	35.5"
Overall Height	44"	44"	44"	44"	43"
Width Between Arms	21"	22"	21"	20"	18"
Floor to Top of Seat	19.5"	20"	19.5"	19"	22"
Seat Depth	21"	21.5"	21"	21"	22"
Seat to Top of Back	28"	28"	28"	29"	27"
Seat to Top of Arm	6.5"	8"	6"	8"	8"
Floor to Back Edge of Seat Raised	27"	28"	27.5"	28.5"	30"
Floor to Front Edge of Seat Raised	23"	24"	23.5"	25"	24"
Seat Angle Raised	23 DEG.	18 DEG.	15 DEG.	15 DEG.	24 DEG.
Floor to Top of Foot Rest Reclined	23"	24"	23.5"	23"	21.5"
Floor to Top of Back Reclined	34.5"	33"	33"	21"	37"
Back Angle Closed	102 DEG.	110 DEG.	115 DEG.	100 DEG.	120 DEG.
Back Angle Reclined	138 DEG.	140 DEG.	145 DEG.	174 DEG.	140 DEG.
Foot Rest Angle Closed	90 DEG.	90 DEG.	90 DEG.	90 DEG.	90 DEG.
Foot Rest Angle Opened	180 DEG.	180 DEG.	180 DEG.	180 DEG.	180 DEG.
Distance Required from Wall Reclined	18"	18"	19"	25"	6"
Extension reclined from Back of Seat to end of Foot Rest	39"	42"	42"	43"	44"
Life Time Warranty, Chair Frame, Lift Frame, Recline Mechanism	YES	YES	YES	YES	1YR
Warranty on Electrical parts	3 YRS.	3 YRS.	3 YRS.	3 YRS.	1 YR.
Pro Rated Years 4 Thru 7	YES	YES	YES	YES	NO
Warranty on Mechanical Labor	3 YRS.	3 YRS.	3 YRS.	3 YRS.	1 YR.
Warranty on Electrical Labor	1 YR.	1 YR.	1 YR.	1 YR.	1 YR.
Removable Back	YES	YES	YES	YES	YES
Customizable *	NO	NO	NO	NO	NO
Left Hand Control Available	YES	YES	YES	YES	YES
Extra Pocket Available	YES	YES	YES	YES	YES
Smart Tek System	YES	YES	YES	NO	NO
Heat Massage Option	YES	YES	YES	YES	YES
Foot Rest Extension Option	NO	NO	YES	YES	NO
MaxiComfort Option	NO	NO	NO	YES	NO



COMFORTER SERIES

	Comforter	Comforter	Comforter	Comforter	Comforter
Model Number	PR-501JP	PR-501S	PR-501M	PR-501L	PR-501T
Size	Jr. Petite	Small	Medium	Large	Tall
Weight Capacity	300LB	300LB	375 LB	375LB	375LB
Recline Positions	3	3	3	3	3
Standard Fabrics	8	8	8	8	8
Fabric Options	36	36	36	36	36
Chaise Seat	YES	YES	YES	YES	YES
Back Type	SEAM	SEAM	SEAM	SEAM	SEAM
Head Rest Cover	YES	YES	YES	YES	YES
Arm Covers	YES	YES	YES	YES	YES
Battery Back Up	YES	YES	YES	YES	YES
Overall Width	30.5"	31.5"	32.5"	34.5"	34.5"
Overall Height	39"	40"	42"	44"	47"
Width Between Arms	19"	20"	21"	23"	23"
Floor to Top of Seat	18"	19"	19.5"	20.5"	22.5"
Seat Depth	18.5"	19.5"	21"	22"	24"
Seat to Top of Back	22.5"	25"	27"	28.5"	30"
Seat to Top of Arm	6.5"	6"	7"	5.5"	6"
Floor to Back Edge of Seat Raised	26"	27.5"	27.5"	29"	31"
Floor to Front Edge of Seat Raised	22"	24"	23.5"	25"	26"
Seat Angle Raised	16 DEG.	15 DEG.	16 DEG.	16 DEG.	18 DEG.
Floor to Top of Foot Rest Reclined	22"	23"	23"	23"	23"
Floor to Top of Back Reclined	32"	32"	32"	33"	36"
Back Angle Closed	115 DEG.	115 DEG.	115 DEG.	115 DEG.	115 DEG.
Back Angle Reclined	145 DEG.	145 DEG.	145 DEG.	145 DEG.	145 DEG.
Foot Rest Angle Closed	90 DEG.	90 DEG.	90 DEG.	90 DEG.	90 DEG.
Foot Rest Angle Opened	180 DEG.	180 DEG.	180 DEG.	180 DEG.	180 DEG.
Distance Required from Wall Reclined	17"	17"	17"	18"	19"
Extension reclined from Back of Seat to end of Foot Rest	38"	39"	42"	43"	45"
Life Time Warranty, Chair Frame, Lift Frame, Recline Mechanism	YES	YES	YES	YES	YES
Warranty on Electrical parts	3 YRS.	3 YRS.	3 YRS.	3 YRS.	3 YRS.
Pro Rated Years 4 Thru 7	YES	YES	YES	YES	YES
Warranty on Mechanical Labor	3 YRS.	3YRS.	3 YRS.	3 YRS.	3 YRS.
Warranty on Electrical Labor	1 YR.	1 YR.	1 YR.	1 YR.	1 YR.
Removable Back	YES	YES	YES	YES	YES
Customizable *	YES	YES	YES	YES	YES
Left Hand Control Available	YES	YES	YES	YES	YES
Extra Pocket Available	YES	YES	YES	YES	YES
Smart Tek System	YES	YES	YES	YES	YES
Heat Massage Option	YES	YES	YES	YES	YES
Foot Rest Extension Option	NO	YES	YES	YES	YES
MaxiComfort Option	NO	YES*	YES*	YES*	YES



COMFORTER WIDES

	Comforter	Comforter	Comforter	Comforter	Comforter
Model Number	PR-501S-XW	PR-501M-26D	PR-501L-26D	PR-501T-28D	PR-502
Size	Wide	Medium	Large	Tall	Wide
Weight Capacity	375 LB	500LB	500LB	500LB	700LB
Recline Positions	3	3	3	3	INFINITE
Standard Fabrics	8	4	4	4	4
Fabric Options	36	36	36	36	36
Chaise Seat	YES	YES	YES	YES	NO
Back Type	SEAM	SEAM	SEAM	SEAM	SEAM
Head Rest Cover	YES	YES	YES	YES	YES
Arm Covers	YES	YES	YES	YES	YES
Battery Back Up	YES	YES	YES	YES	YES
Overall Width	34.5"	38.5"	38.5"	41"	45"
Overall Height	40"	42"	44"	48"	47"
Width Between Arms	23"	27.5"	27.5"	29.5"	33"
Floor to Top of Seat	18"	19.5"	20.5"	22.5"	21"
Seat Depth	19.5"	21.5"	22.5"	24"	22"
Seat to Top of Back	25"	27"	28.5"	30"	29"
Seat to Top of Arm	6"	7"	5.5"	6"	7"
Floor to Back Edge of Seat Raised	27.5"	27.5"	29"	31"	31"
Floor to Front Edge of Seat Raised	24"	23.5"	25"	26"	27"
Seat Angle Raised	15 DEG.	16 DEG.	16 DEG.	18 DEG.	12 DEG.
Floor to Top of Foot Rest Reclined	23"	23"	23"	23"	21"
Floor to Top of Back Reclined	32"	32"	33"	36"	39"
Back Angle Closed	115 DEG.	115 DEG.	115 DEG.	115 DEG.	110 DEG.
Back Angle Reclined	145 DEG.	145 DEG.	145 DEG.	145 DEG.	150 DEG.
Foot Rest Angle Closed	90 DEG.	90 DEG.	90 DEG.	90 DEG.	90 DEG.
Foot Rest Angle Opened	180 DEG.	180 DEG.	180 DEG.	180 DEG.	180 DEG.
Distance Required from Wall Reclined	17"	17"	18"	19"	20"
Extension reclined from Back of Seat to end of Foot Rest	39"	42"	43"	45"	45"
Life Time Warranty, Chair Frame, Lift Frame, Recline Mechanism	YES	YES	YES	YES	2 YRS.
Warranty on Electrical parts	3 YRS.	3 YRS.	3 YRS.	3 YRS.	1 YR.
Pro Rated Years 4 Thru 7	YES	YES	YES	YES	NO
Warranty on Mechanical Labor	3 YRS.	3 YRS.	3 YRS.	3 YRS.	2 YRS.
Warranty on Electrical Labor	1 YR.	1 YR.	1 YR.	1 YR.	1 YR.
Removable Back	YES	YES	YES	YES	YES
Customizable *	YES	YES	YES	YES	YES
Left Hand Control Available	YES	YES	YES	YES	YES
Extra Pocket Available	YES	YES	YES	YES	YES
Smart Tek System	YES	YES	YES	YES	YES (1/05)
Heat Massage Option	YES	YES	YES	YES	YES
Foot Rest Extension Option	YES	YES	YES	YES	NO
MaxiComfort Option	NO	NO	NO	NO	NO



VALUE SERIES

	EZ Seat	Capri	Slimline	Monarch	Monarch
Model Number	PR-150	PR-CAPRI	PR-100	PR-355M	PR-355L
Size	Small	Medium	Medium	Medium	Large
Weight Capacity	250LB	325LB	325LB	375LB	375LB
Recline Positions	NA	2	2	3	3
Standard Fabrics	2	7	3	7	7
Fabric Options	0	7	0	36	36
Chaise Seat	NO	NO	NO	NO	NO
Back Type	STRAIGHT	BUTTON	SEAM	SEAM	SEAM
Head Rest Cover	NO	YES	YES	YES	YES
Arm Covers	NO	YES	YES	YES	YES
Battery Back Up	YES	YES	YES	YES	YES
Overall Width	24"	31"	29"	32"	34"
Overall Height	41.5"	41"	41"	42"	44"
Width Between Arms	19.25"	20"	20.5"	21"	22"
Floor to Top of Seat	18.5"	18"	19"	18.5"	18.5"
Seat Depth	18.5"	19"	19"	21.5"	22"
Seat to Top of Back	24.5"	27"	27"	28"	30"
Seat to Top of Arm	8.5"	7"	6.5"	7"	7"
Floor to Back Edge of Seat Raised	25.5"	27"	26"	27"	27"
Floor to Front Edge of Seat Raised	25.5"	22"	22.5"	23"	23"
Seat Angle Raised	180 DEG.	16 DEG.	18 DEG.	16 DEG.	16 DEG.
Floor to Top of Foot Rest Reclined	NA	23"	23"	23"	22"
Floor to Top of Back Reclined	NA	35"	36"	33"	33"
Back Angle Closed	NA	110 DEG.	100 DEG.	110 DEG.	110 DEG.
Back Angle Reclined	NA	135 DEG.	125 DEG.	140 DEG.	140 DEG.
Foot Rest Angle Closed	NA	90 DEG.	90 DEG.	90 DEG.	90 DEG.
Foot Rest Angle Opened	NA	160 DEG.	160 DEG.	176 DEG.	176 DEG.
Distance Required from Wall Reclined	0"	11"	11"	17"	18"
Extension reclined from Back of Seat to end of Foot Rest	NA	40"	38"	43"	43"
Life Time Warranty, Chair Frame, Lift Frame, Recline Mechanism	1 YR.	YES	YES	YES	YES
Warranty on Electrical parts	1 YR.	1 YR.	1 YR.	2 YRS.	2 YRS.
Pro Rated Years 4 Thru 7	NO	NO	NO	NO	NO