

S I G N A T U R E S E R I E S



# Rise



## Space Saver

Plush, overstuffed support combined with our highly engineered space saving recline system requiring only 6" of space for a full recline.



# Space Saver

S I G N A T U R E S E R I E S

Model Number	PR-906
Size	Medium
Weight Capacity	250 LB
Recline Positions	3
Standard Fabrics	3
Fabric Options	3
Chaise Seat	YES
Back Type	SEAM
Head Rest Cover	NO
Arm Covers	NO
Battery Back Up	YES
Overall Width	35.5"
Overall Height	43"
Width Between Arms	18"
Floor to Top of Seat	22"
Seat Depth	22"
Seat to Top of Back	27"
Seat to Top of Arm	8"
Floor to Back Edge of Seat Raised	30"
Floor to Front Edge of Seat Raised	24"
Seat Angle Raised	24°
Floor to Top of Foot Rest Reclined	21.5"
Floor to Top of Back Reclined	37"
Back Angle Closed	120°
Back Angle Reclined	140°
Foot Rest Angle Closed	90°
Foot Rest Angle Opened	180°
Distance Required from Wall Reclined	6"
Extension Reclined from Back of Seat to End of Foot Rest	44"
Life Time Warranty, Chair Frame, Lift Frame, Recline Mechanism	1 YR.
Warranty on Electrical Parts	1 YR.
Pro Rated Years 4 Thru 7	NO
Warranty on Mechanical Labor	1 YR.
Warranty on Electrical Labor	1 YR.
Removable Back	YES
Customizable	NO
Left Hand Control Available	YES
Extra Pocket Available	YES
Smart Tek System	NO

Sizes may vary depending on fabric, filling material, upholstery and carpet thickness. All measurements taken on concrete floor using levels and metal rulers.

Hand-blown, Overstuffed, Conjugated Microfiber Dacron Fill  
**SUPER SUPPORTIVE**

Northern Maple Wood  
• Double-doweled at Stress Points  
• Reinforced Miter Blocks  
• Stapled and Glued  
**STRONGEST IN THE INDUSTRY**

2" of Block Foam  
**EXTRA COMFORT**

Quick Release Hand Control  
**EASY DISCONNECT**


Rear Landing Legs  
**EXTRA STABILITY**

5 Heavy Gauge Springs  
**LONG LASTING FOR INCREDIBLE COMFORT**

Widest Fabric Selection in Industry  
**FITS IN ANY HOME DÉCOR**

10 Gauge Steel Scissor Mechanism  
**STRONG & LONG LASTING**

Adjustable Anti-Skid/Scrape Feet  
**STABILITY & SUPPORT**



501 M shown



501 M shown



501 M shown



501 M shown

## Golden Craftsmanship

Each luxury lift & recline chair is hand crafted in our facility using decades of motion furniture experience.

## Added Comfort and Safety

Designed to perform year after year using advanced technology backed by the best warranty.

## Smart Engineering

Smooth, quiet lifting system offers a secure lift for the best performance.



Battery Backup System

## Safety System

Batteries not included. Safety System is designed to bring chair from recline to seated position in the event of a power outage.



Space Saver shown in Atlantic fabric



Literature is current at time of printing. Company reserves the right to make changes.



Our chairs meet some of the toughest safety standards ever. Our entire chair meets ETL (Electrical Testing Lab) standards to ensure electrical safety. ETL is a nationally recognized, third party testing facility.



**GOLDEN**  
technologies