

APEXLIFT

Heavy Duty Electric Patient Lift



BENEFITS

- 24 VDC Battery Operated
- Safe Working Load 600lbs
- Extra wide Spreader Bar for Patient Comfort
- Adjustable Spreader Bar Angle for Positioning
- Low Battery/Overload Safety Warning Devices

TECHNICAL INFORMATION

- Safe Load Capacity : 600 lbs
- Hook width: 27.5" (700mm)
- Over Length : 50" (1270 mm)
- Unit Weight : 137 lbs(62kgs)
- Front Caster Size: 4" (102mm)
- Rear Braked Caster Size: 4" (102mm)
- Spreader Bar Min Height : 27" (690mm)
- Spreader Bar Max Height : 71" (1810mm)
- Widest Internal Base Width: 39" (660mm)
- Narrowest External Base Width : 26" (990mm)
- Volt: 24VDC



Apex Dynamics Healthcare Products, L.L.C.