


# EZ POOL LIFT™



## SPECIFICATIONS

THE AQUA CREEK SURE HANDS LIFT SHALL BE A HAND POWERED SPA OR IN-GROUND POOL LIFT. THE LIFT SHALL BE STAINLESS STEEL CONSTRUCTION WITH A LIFTING CAPACITY OF 400 LBS. THE LIFT SHALL HAVE AN OPERATING RANGE OF 53" AND ADJUSTABLE TO A RANGE OF 70". THE SURE HANDS LIFT CAN BE SETBACK UP TO 20" FROM THE POOL EDGE WITH A 360° MANUAL ROTATION. THE UNIT SHALL INCLUDE A 6" DEEP ANCHOR WITH AN ANTI-TAMPER LID. THE UNIT IS PORTABLE AND CAN BE REMOVED WITHOUT TOOLS. A 5 YEAR STRUCTURAL AND 1 YEAR SLING WARRANTY IS STANDARD.

MAST: 1.9" O.D. .145" WALL STAINLESS STEEL TUBE TYPE 304 WITH REINFORCING FIN  
 BOOM: 1.9" O.D. .145" WALL STAINLESS STEEL TUBE TYPE 304 WITH REINFORCING FIN  
 RAM: SELF-CONTAINED POWER-PACKER HYDRAULIC RAM WITH NITROTEC FINISH  
 SLING: 100% POLYESTER MESH RATED FOR PATIENT LIFTING DEVICES  
 CHAIN: STAINLESS STEEL TYPE 304  
 FINISH: UV RESISTANT POLYESTER POWDER-COATING

PRODUCT:	<b>F-003PL</b>	
DWG #:	<b>F-003PL</b>	
SCALE:	<b>1:18</b>	DATE: <b>07/29/05</b>
DRAWN BY:	<b>DSE</b>	REVISION: <b>A</b>