


POWER - EZ POOL LIFT™


THE EZ POOL LIFT WILL CLEAR A 26" ABOVE GROUND POOL OR SPA.



SPECIFICATIONS

THE AQUA CREEK EZ POOL LIFT-BATTERY SHALL BE A BATTERY POWERED SPA OR IN-GROUND POOL LIFT. THE LIFT SHALL BE STAINLESS STEEL CONSTRUCTION WITH A LIFTING CAPACITY OF 400 LBS. THE LIFT SHALL HAVE AN OPERATING RANGE OF 53" AND ADJUSTABLE TO A RANGE OF 70". THE EZ POOL -BATTERY CAN BE SETBACK UP TO 20" FROM THE POOL EDGE WITH A 360° MANUAL ROTATION. THE UNIT SHALL INCLUDE A 6" DEEP ANCHOR WITH AN ANTI-TAMPER LID. THE UNIT IS PORTABLE AND CAN BE REMOVED WITHOUT TOOLS. A 5 YEAR STRUCTURAL AND 1 YEAR SLING WARRANTY IS STANDARD.

MAST: 1.9" O.D. .145" WALL STAINLESS STEEL TUBE TYPE 304 WITH REINFORCING FIN
 BOOM: 1.9" O.D. .145" WALL STAINLESS STEEL TUBE TYPE 304
 RAM: LINAK LA34 LINEAR ACTUATOR
 CONTROLLER: WATERPROOF LINAK CONTROLLER
 BATTERY: 24 DC RECHARGABLE BATTERY PACK
 SLING: 100% POLYESTER MESH RATED FOR PATIENT LIFTING DEVICES
 CHAIN: STAINLESS STEEL TYPE 304
 FINISH: UV RESISTANT POLYESTER POWDER-COATING



PRODUCT:	F-003PEZ	
DWG #:	F-003PEZ	
SCALE:	1:18	DATE: 07/26/05
DRAWN BY:	MBT	REVISION: A

TRITON POOL LIFT

BY WORLD WIDE SEATING



The Triton Pool Lift is has a weight capacity up to 400 pounds and is easily portable for storage and/or relocation to another pool area. Quality stainless steel construction and a UV resistant powder coating finish ensure long life and corrosion resistance. Anchor included. Ships UPS 50 lbs. The Triton Lift comes with a 5 year Warranty

#TRI7HYDR-400 Triton Hydraulic Power Lift

#TRI7ELECT-400 Triton Electric Power

Includes Sling Small, Medium, Large or Extr-Large

Includes Concreat Anchor

Optional Surface Mount Anchor or Wood Deck Anchor