


RMT AQUATICS



NOTES:

1. INSTALLATION TO BE COMPLETED IN ACCORDANCE WITH MANUFACTURER'S SPECIFICATIONS.
 2. DO NOT SCALE DRAWINGS.
 3. LIFTING CAPACITY 300 LB.
 4. TOTAL TRAVEL DISTANCE WITH ACTUATOR IS 58".
 5. SEAT DEPTH 18"-20" BELOW WATER LINE.
 6. 359° TOTAL ROTATION.
 7. CONTRACTORS NOTE: FOR PRODUCT AND PURCHASING INFORMATION VISIT www.CADdetails.com/info
- REFERENCE NUMBER 577-006



SPLASH! EXTENDED REACH SEMI-PORTABLE POOL LIFT

COMPONENT IDENTIFICATION, SIDE VIEW

577-006

PROTECTED BY COPYRIGHT - 05/27/05

www.CADdetails.com