

# ErgoStand



# Ergostand

*For patients with some strength, transfers will be fast and easy with the Ergostand.*



*Your patient is comfortably secure in our band sling with a velcro waist enclosure reinforced with 2 belts.*



*Our handset allows you to work beside your patient and have an extra set of controls on the control box.*



*Superior back support and anti-slip material ensure your patient is comfortable and safe.*



*The fastest and most efficient way to dress, clean and transfer patients who need additional support.*



*Your patient will feel comfortable and safe with the padded knee support and anti-slip foot base.*

*Even though the Ergostand is extra sturdy, it's still small enough to easily fit through doorways and tight spaces.*



*Obstacles are no longer a worry with the motorized base.*




# Ergostand

## FEATURES

Lifting capacity	180 kg (400 lbs.)
Soft start & stop	YES
Motorized base	YES
2-speed vertical movement	YES
Multipositioned handles close to caregiver	YES
Directional wheel guide	YES
Emergency stopping device	YES
Padded handles and leg support	YES
CSA 601.1, UL 2601.1, EN 60601.1	YES
Integrated charging system	YES
Visual and audible low battery indicator	YES


## ERGOFIT SLINGS

All slings are designed to be installed and removed quickly.



### HYGIENIC SLING

Ideal for transferring to the toilet. It also permits an easy way to change diapers.



### BAND SLING

Perfect for the Ergostand. Velcro waist with 2 additional belts, back supports and anti-slip material for the back.

## SLING SIZES AND PRODUCT NUMBERS

Size	SMALL	MEDIUM	LARGE
<b>Suggested weight:</b>	20 - 46 kg (45 - 100 lbs.)	46 - 95 kg (100 - 210 lbs.)	95 - 250 kg (210 - 550 lbs.)
<b>Suggested height:</b>	120 - 150 cm (4' - 4' 11")	151 - 180 cm (5' - 5' 11")	181 cm + (6' +)
<b>Hygienic</b>	THY-S	THY-M	THY-L
<b>Band</b>	TST-S	TST-M	TST-L

For a demonstration or more information:  
1-800-868-0441

2001 Tanguay Street  
Magog, Quebec, Canada J1X 5Y5  
Tel.: (819) 868-0441  
Fax: (819) 868-2249  
Internet: www.bhm-medical.com  
E-mail: bhm@bhm-medical.com