

Pride
Mobility Products Corp.

Silver Star Lifts ★
A Division of Pride Mobility Products Corp.



BACKPACKER INTERIOR LIFT

Your on-the-go lifestyle demands that you bring your scooter or power chair with you. Now, Pride Mobility Products® Corp., the number one power chair and scooter manufacturer, offers the lift you need to take your mobility product on the road. Backed by the service you have come to expect and trust from Pride®, with the full support of our expert technical service staff, a Silver Star™ Lift is the natural complement to your Pride Scooter or Power Chair.

Take your mobility on the road with the Silver Star Backpacker interior lift from Pride. It's the ultimate solution for lifting your power chair or scooter into your full size van or minivan for transport. Like all Silver Star Lifts from Pride, the Backpacker features all-steel construction and a protective finish, along with an industrial rated 12 volt motor and optional battery pack, and easy-to-use controls to combine rugged durability with safe, simple and reliable operation.

Silver Star™ Lifts ★

A Division of **Pride** Mobility Products Corp.

BACKPACKER INTERIOR LIFT



Backpacker in lowest position.



Backpacker raising.



Backpacker loaded.

SPECIFICATIONS

- Base weight lift** 106lbs. ¹
- Scooter / power chair overall weight** 325 lbs.
- Maximum product width** 26"
- Maximum product wheelbase** 35"
- Maximum travel distance from the ground** 26"
- Motor** Sealed 12 volt DC
- Warranty** 3-year transferable
- Optional battery pack** 12 volt, on-board
- Optional charger** 110 volt



Can easily lift the following Jazzy® Power Chairs ²

- Jazzy® 600
- Jazzy® 610
- Jazzy® 1103 Ultra
- Jazzy® 1107
- Jazzy® 1113 ATS
- Jazzy® 1121
- Jazzy® 1122
- Jazzy® 1133
- Jazzy® 1143 Ultra

Can easily lift the following Quantum® Power Chairs ²

- Quantum® 600
- Quantum® 610
- Quantum® 1107
- Quantum® 1121
- Quantum® 1402

Can easily lift the following Pride® Scooters ²

- Celebrity® X
- Go-Go® Ultra
- Go-Go Elite Traveller™
- Legend®
- Maxima
- Rally®
- Revo®
- Sonic®
- Sundancer®
- Victory®

Can easily lift the following Jet™ Power Chairs ²

- Jet™ 2
- Jet™ 2HD
- Jet™ 3 Ultra
- Jet™ 7
- Jet™ 10
- Jet™ 10 Ultra
- Go-Chair™

¹ Weight does not include any accessories
² Vehicle dependent

Note: The Backpacker can lift all of the listed power chairs with contoured, Synergy® or Specialty seating. Please call Customer Service when requesting information about lifting with a power positioning system.

Silver Star™ Lifts
 A Division of **Pride** Mobility Products Corp.