



Revo®

3 & 4-wheel scooters



The Revo® delivers “full-sized” scooter style, performance, and features in a compact, lightweight design that easily disassembles for transport while remaining highly affordable.

Introducing the Revo 4-wheel.

The Revo provides sleek styling, super portability, and now the enhanced stability of a 4-wheel base.



Pride
Mobility Products Corp.

Specifications

Model number(s)	3-wheel: SC63 4-wheel: SC64
Weight capacity	300 lbs.
Maximum speed	Up to 4.5 mph
Per-charge range	Up to 20 miles
Ground clearance	2.25"
Turning radius	3-wheel: 38" 4-wheel: 50"
Overall length	3-wheel: 43" 4-wheel: 44.75"
Overall width	3-wheel: 21.5" 4-wheel: 21.5"

Tires

Front	2" x 8" solid
Rear	3" x 9" solid
Anti-tips	Rear
Standard seating	Compact foldable molded plastic
Power elevating seat	NA
Suspension	NA
Drive train	Sealed transaxle, 24 volt DC motor
Dual braking system	Electronic, regenerative and electromechanical
Battery charger	Off-board, 5 amp
Battery requirements	(2) 12V deep-cycle (AGM or gel cell)
Battery size(s)	(U-1)

Weight

Battery	24.5 lbs. each
Total w/o battery packs	3-wheel: 106.5 lbs. 4-wheel: 125 lbs.
Heaviest piece when disassembled	3-wheel: 48 lbs. (rear) 4-wheel: 51.5 lbs. (front)

Colors



Blue



Red

The information contained herein is correct at the time of publication; we reserve the right to alter specifications without prior notice. Speed & range vary with user weight, terrain type, battery charge, battery condition and tire pressure.

Benefits

- Innovative interlocking drop-in battery packs for easy removal/installation
- Easily disassembles into 5 lightweight, easy-to-manage pieces
- 300 lbs. weight capacity
- Controls at your fingertips
- Easy-to-use external freewheel lever
- State-of-the-art, microprocessor-based motor controller
- Seat offers height adjustment, adjustable 360° seat rotation, and flip-up armrests for optimal user comfort
- Adjustable tiller for maximum comfort
- Headlight and taillight for added safety
- On-board accessory bracket
- Easily accessible charger port and circuit breaker
- 3-year limited warranty

Options

- Forearm crutch holder
- Lap belt
- Cane/crutch holder
- Wishbone crutch holder
- Oxygen tank holder
- Cup holder
- Rear basket
- Safety flag
- Walker holder
- Heavy-duty dust cover
- Saddlebag
- Quad cane holder
- Cane, clip attachment

