



888B-3
Sunrunner-3



888-3L
Sunrunner-3 Deluxe



888B-4
Sunrunner-4



888-4L
Sunrunner-4 Deluxe

PRICE LIST & ORDER FORM Effective July 2004

For easy ordering and information
contact Customer Service Toll Free at:

800-743-0772 Fax: 877-797-7081

Name _____ Phone No _____
 Acct No _____ PO # _____ Date _____
 Address _____ City/State _____ Zip _____
 Ship To _____ Phone No _____
 Address _____ City/State _____ Zip _____

Stock For your convenience, order by stock number (e.g. BAT-U1)

BASE SCOOTER

<input type="checkbox"/>	888B-3 Sunrunner 3	Burgundy/Blue	E1230	\$2,195.00
<input type="checkbox"/>	888B-4 Sunrunner 4	Burgundy/Blue	E1230	\$2,395.00
<input type="checkbox"/>	888-3L SUNRUNNER 3 DELUXE	Burgundy	E1230	\$2,295.00
<input type="checkbox"/>	888-4L SUNRUNNER 4 DELUXE	Burgundy	E1230	\$2,495.00

BATTERY

<input type="checkbox"/>	BAT-U1	12V 35AMP	x2	Standard
--------------------------	--------	-----------	----	----------

SEATS (W"xD"xH")

<input type="checkbox"/>	DXV-171/2	17½" Deluxe Captain Seat Vinyl-888B-3 & B-4	Standard
<input type="checkbox"/>	104018-00102	18"x18" Pillow-Top Recline Seat Vinyl-888-3L & 4L	Standard
<input type="checkbox"/>	DXC-171/2	Upgrade 17½" Deluxe Captain Seat Grey Cloth	\$100.00
<input type="checkbox"/>	DXV-18*20R	Upgrade 18"x20" Recline Deluxe Cpt Seat Vinyl	\$100.00
<input type="checkbox"/>	DXV-20*22	Upgrade 20"x22" Recline Capt Grey Vinyl	\$100.00

ARM REST ASSY

<input type="checkbox"/>	104300-88800	(L) Arm Rest Assy w/ Arm Pad	Standard
<input type="checkbox"/>	104400-88800	(R) Arm Rest Assy w/ Arm Pad	Standard

TIRES

<input type="checkbox"/>	TAR26085	Rear	10" Pnuematic Tires-all Models	x2	Standard
<input type="checkbox"/>	205102-88100	Rear	10" Tire Tubes-all Models	x2	Standard
<input type="checkbox"/>	TAF26085	Front	10" Pnuematic Tires-B-4 & 4L	x2	Standard
<input type="checkbox"/>	205102-88100	Front	10" Tire Tubes-B-4 & 4L	x2	Standard
<input type="checkbox"/>	TAF280250	Front	9" Pnuematic Tires-B-3 & 3L	x1/x2	Standard
<input type="checkbox"/>	205202-77100	Front	9" Tire Tubes-B-3 & 3L	x2	Standard
<input type="checkbox"/>	STU8-9(K)	F&R	888B-3 Solid Tire Upgrade		\$200.00
<input type="checkbox"/>	STU8WS(K)	F&R	888-3L Solid Tire Upgrade		\$280.00
<input type="checkbox"/>	STU8A(K)	F&R	888B-4/888-4L Solid Tire Upgrade		\$280.00

CHARGER

<input type="checkbox"/>	N/A	4 amp Charger On-Board (888-3L & 888-4L)	Standard
<input type="checkbox"/>	109201-88306	3 amp Charger Off-Board (888B-3 & 888B-4)	Standard

ACCESSORIES

<input type="checkbox"/>	109901-8810*	Front Wicker Basket	Standard
<input type="checkbox"/>	102304-67200	Crutch Holder Bracket	\$40.00
<input type="checkbox"/>	CCH-003	Cane Holder	\$150.00
<input type="checkbox"/>	CH-02	Oxygen Tank Holder (Metal)	\$150.00
<input type="checkbox"/>	WH-0185	Walker Holder	\$150.00
<input type="checkbox"/>	OXY-03	Mesh Oxygen Tank Holder	\$60.00
<input type="checkbox"/>	RB-0125	Rear Wire Mesh Basket	\$170.00
<input type="checkbox"/>	300909-03	Safetv Belt	\$34.00

All specifications are subject to change without notice.
 Shoprider Mobility Products, Inc® reserves the rights of any changes to the unit:

- 1) Include the anti-tip wheel or rear castor.
- 2) Driver Weight may exceed weight of the unit; speed must be reduced when turning.
- 3) The actual driving range varies with the factors shown below:
 - a) The weight of user.
 - b) Ground surface and terrain.
 - c) Battery capacity and conditions.
 - d) Type of charger.
 - e) Ambient Temperature.
 - f) The way and/or method of driving.
 - g) If the battery and mechanical moving parts are fully broken in.

TOTAL \$

MIDSIZE SPECIFICATIONS

SUNRUNNER-3 (888B-3)

Base LxWxH*	45"x24"x38"
*Ground to Seat Top	
Ground Clearance	3"
User Weight Capacity	300 Lbs
Total Unit Weight	165 Lbs
Battery Module Weight (U1)	25 Lbs x2
Number, Size Tire (front/rear)	1,9"/2,10"
BAT-U1 Module Capacity	12v35Ah x2
Charger-Off Board	3 amp
Turning Radius	37"
Range	up to 25 miles
Maximum Speed	5 mph
Maximum Climbing Angle	14% / 8°

SUNRUNNER-4 (888B-4)

Base LxWxH*	49"x22"x38"
*Ground to Seat Top	
Ground Clearance	3"
User Weight Capacity	300 Lbs
Total Unit Weight	180 Lbs
Battery Module Weight (U1)	25 Lbs x2
Number, Size Tire (front/rear)	2,10"/2,10"
BAT-U1 Module Capacity	12v35Ah x2
Charger-Off Board	3 amp
Turning Radius	50"
Range	up to 25 miles
Maximum Speed	5 mph
Maximum Climbing Angle	14% / 8°

SUNRUNNER-3 Deluxe(888-3L)

Base LxWxH*	48-51"x24"x41"
*Ground to Seat Top	
Ground Clearance	3"
User Weight Capacity	300 Lbs
Total Unit Weight	160 Lbs
Battery Module Weight (U1)	25 Lbs x2
Number, Size Tire (front/rear)	2,9"/2,10"
BAT-U1 Module Capacity	12v35Ah x2
Charger-On Board	4 amp
Turning Radius	39"
Range	up to 24 miles
Maximum Speed	5 mph
Maximum Climbing Angle	14% / 8°

SUNRUNNER-4 Deluxe(888-4L)

Base LxWxH*	49-52"x24"x41"
*Ground to Seat Top	
Ground Clearance	3"
User Weight Capacity	300 Lbs
Total Unit Weight	176 Lbs
Battery Module Weight (U1)	25 Lbs x2
Number, Size Tire (front/rear)	2,10"/2,10"
BAT-U1 Module Capacity	12v35Ah x2
Charger-On Board	4 amp
Turning Radius	53"
Range	up to 24 miles
Maximum Speed	5 mph