

GRANSTAND® II ORDER FORM

Client Name:			PO#	Date:
Age:	Height	Weight	Dealer:	
Address:			Dealer Contact:	
City, St, Zip:			Address:	
Phone:			City, St, Zip:	
Facility:			Phone:	
Therapist:		Phone	Fax:	

GRANSTAND® II: Hydraulic Assisted Standing Frame

All Granstand Standing Systems include our **NO TRANSFER DESIGN**, eliminating the risk of injury from transfers and lifts and allowing the client to be lifted directly from their mobility product. The Granstand II's innovative base design allows lifting from the wheelchair, a bed, a mat table or other surfaces and re-seating them safely and properly onto any of these surfaces

GRANSTAND® II Frame Styles: Standard Features include: Expandable multi-functional transport base with 3" plastic casters, depth adjustable chest pad, adjustable angle clear table top, width and depth adjustable foot positioners, pull-to-place seat sling and depth adjustable knee pad.

- Granstand II fits client height: 5'6" – 6'4". (Standard)
 - Granstand II fits client height: 4'6" – 6'5" (requires one of the mobile foot systems listed below)
- Pump Handle** Right Left **Frame Color Choices:** Black Red Blue Granite

Sling Style: All slings include our easy **pull-to-place** strap system.

- Web Sling (Standard):** Easy to use **pull-to-place** sling in 2" black nylon webbing.
- Sheepskin Sling Cover:** Genuine sheepskin cover with Velcro attachments.
- Pocket Seat Sling Cover:** Open pocket design allows for additional padding of the sling as desired.
- Vinyl Seat Sling:** Vinyl easy to use **pull-to-place** seat sling allows greater surface contact and ease of cleaning.



Lower Extremity Positioning:

- Standard Foot System:** Allows for width and depth adjustment. (*Granstand II is not mobile in use with this item*)
- Height-Adjustable Mobile Foot Support System:** Adjusts in width, depth and height allowing for proper adjustment of the foot positioners to the height of the footplates of the wheelchair. Supplies maximum growth potential and height adjustment plus allows the Granstand II to be moved while in use.
- Multi-Adjustable Mobile Foot Support System:** Adjusts in width, depth, height and angle for precise foot position, plus allows for individual adjustments for planter-flexion or dorsi-flexion contractures. Supplies maximum growth potential and height adjustment plus allows the Granstand II to be moved while in use.
- Standard Knee System:** Depth and height adjustable one piece padded knee support
- Multi-Adjustable Knee System:** Quick adjust mechanism easily adjusts in height, width, and depth for optimum lower body positioning. Independent width positioning allows for proper placement of the knee pad and fast adjustment for multiple client use.
- Roho Knee Pads:** Used with the Multi-Adjustable Knee System this option offers additional padding for clients with special positioning needs such as contractures.



(Not pictured)



PRIME ENGINEERING 4838 W JACQUELYN AVE, #105, FRESNO, CA 93722 1-800-827-8263

A Division of Axiom Industries Inc. Fax (800) 800-3355 E-mail info@Primeengineering.com Web site: www.primeengineering.com

